

Hartmann's Aural Speculum



This instrument is used to straighten the external auditory canal to examine the tympanic membrane.

Diagnostic uses.

This instrument is used to examine the external auditory canal and tympanic membrane.

Therapeutic uses of Hartmann's Aural Speculum.

Hartmann's Aural Speculum has the following uses for therapeutic purposes:

1. Aural toileting in case of CSOM (Chronic Suppurative Otitis Media).
2. Removal of foreign body.
3. Wax removal.
4. During surgeries like: (a) Myringotomy, (b) Myringoplasty and (d) Stapedectomy.
5. Can be used as dilator in stenosis of external auditory canal.

Know it in more detail with FAQs for Viva at:
[Dr. Factopedia](#)

ENT Instrument by Dr. Factopedia

